

BINAURALS FO MAGISTER II + TUBING BLACK (WITHOUT TIPS)




Ref. 506982



 30/08/2024 18:16:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BINAURALS FO MAGISTER II + TUBING BLACK (WITHOUT TIPS)



Ref. 506982

Description

Tubing with matte black lyre for SPENGLER MAGISTER II stethoscope. (without tips)

Features

CE Marking	NONE	Medical device	No
Bar code	3700446034337	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.06
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:16:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr