

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507308



CE



30/08/2024 18:24:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507308

=====

Description

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass. Elegant and streamlined for easy handling. High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips. Ideal for monitoring blood pressure and general auscultation.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446026523	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.393
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Red	Packed unit / Quantity	1

30/08/2024 18:24:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr