



SPENGLER STETHOSCOPE PULSE II SINGLE, PAEDIATRIC



Ref. 507704



CE



30/08/2024 18:37:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPENGLER STETHOSCOPE PULSE II SINGLE, PAEDIATRIC



Ref. 507704

Stethoscope, binaural, complete: A mechanical listening device designed for listening to sounds from the heart and lungs from paediatric patients.

Description

HIGH ACOUSTIC QUALITY

- Aluminium single head chestpiece
- High performance single tubing
- Aluminium binaural
- High performance diaphragm

COMFORT

- Non-chill ring
- Silicone eartips
- New binaural design better adapted to human anatomy

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029906	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	3
Packed unit / Length (cm)	24	Packed unit / Weight	0.215
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:37:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr