

NON CHILL RING (RED) FOR BELL OF DUAL PULSE II SPENGLER



Ref. 507945

 30/08/2024 18:50:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

NON CHILL RING (RED) FOR BELL OF DUAL PULSE II SPENGLER



Ref. 507945

Description

Spare part for SPENGLER DUAL PULSE® II stethoscopes.
Coral-colored anti-cold ring for small bell.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030094	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:50:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr