

## NON CHILL RING (BLACK) FOR CARDIOPRESTIGE II STETHOSCOPE SPENGLER




Ref. 507972



 30/08/2024 18:53:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## NON CHILL RING (BLACK) FOR CARDIOPRESTIGE II STETHOSCOPE SPENGLER



Ref. **507972**

### Description

Spare part for SPENGLER Cardio Prestige® II stethoscopes.  
Anti-cold ring for adult roof in black.

### Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | NONE          | Medical device            | No       |
| Bar code                           | 3700446030315 | HS Code                   | 90181990 |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.001    |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | NA            | Max temperature           | NA       |
| Minimum humidity                   | NA            | Maximum humidity          | NA       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Packed unit / Quantity             | 1             |                           |          |

 30/08/2024 18:53:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)