

FLOATING DIAPHRAGM (SMALL) FOR MAGISTER II DUAL SPENGLER



Ref. 507975S



 30/08/2024 18:54:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FLOATING DIAPHRAGM (SMALL) FOR MAGISTER II DUAL SPENGLER



Ref. 507975S

Description

Spare part for stethoscopes MAGISTER® II from SPENGLER.
Multifrequency fiberglass membrane for pediatric pavilion.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446033187	HS Code	90181990
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.002
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:54:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr