

FLOATING DIAPHRAGM FOR CARDIOPRESTIGE II SPENGLER




Ref. 507976



 30/08/2024 18:55:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FLOATING DIAPHRAGM FOR CARDIOPRESTIGE II SPENGLER



Ref. 507976

Description

Spare part for SPENGLER Cardio Prestige® II stethoscopes.
Multi-frequency fiberglass membrane for adult pavilion.

Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3700446030353 | HS Code | 90181990 |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.001 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

30/08/2024 18:55:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr