

NON CHILL RING (GREEN) FOR PULSE II PEDIATRIC STETHOSCOPE SPENGLER




Ref. 507983

 30/08/2024 18:56:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

NON CHILL RING (GREEN) FOR PULSE II PEDIATRIC STETHOSCOPE SPENGLER



Ref. 507983

Description

Spare part for SPENGLER pediatric PULSE® II stethoscopes.
Anti-cold ring for pediatric pavilion.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030421	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 18:56:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr