

LIFE LINK TUBE FOR MOBI

Ref. 513720



30/08/2024 19:13:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LIFE LINK TUBE FOR MOBI



Ref. 513720

Description

Spare part: LIFE LINK tubing for SPENGLER MOBI blood pressure monitor.
Patented tubing, magnetized for a cuff change in an instant.

Features

CE Marking	NONE	Medical device	Class Im (measuring function)
Bar code	3700446028046	HS Code	90189010
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.199
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:13:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr