



## AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER




Ref. SSDEFI03



 30/08/2024 23:55:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER



Ref. SSDEFI03

**Simple automatic Automated External Defibrillator (AED) within everyone's reach**

## Description

The FRED® PA-1 is perfectly suited for use on both adult and child patients thanks to its dedicated electrodes.

Suitable for use in healthcare institutions, ERPs: it is by design a device intended for the widest audience.

- Electric shock delivery light indicator (for fully automatic version)
- Simple and uncluttered user interface, speech synthesis in French- Has various cardiac massage aid functions to guide the user to perform quality CPR at the right frequency
- Instantly turned on as soon as the hood is opened- Equipped by default with a metronome allowing the user to know the appropriate rhythm of cardiac massage
- Performs periodic self-tests on a daily or weekly basis
- By default in French (Trilingual option possible)
- Can record up to 2 hours of data and 500 events
- Pair of pre-connected defibrillation electrodes housed in the hood of the AED with the first aid kit\* (location for a second pair of electrodes under the AED)
- Indicator lights visible at more than 15m
- ABS V0 box (flame retardant) medical grade- Maximum noise level in operation: 85 db
- Comes standard with 1 battery: LiMnO2 15V, 2.8Ah

\*First aid kit: 1 pair of emergency scissors, 1 disposable razor, 1 pair of vinyl gloves, 1 ventilation mask, 2 sterile compresses, 1 towel to wipe the victim

 30/08/2024 23:55:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## AUTOMATIC DEFIBRILLATOR FRED PA-1 SCHILLER



Ref. **SSDEFI03**

### Features

CE Marking	CE0459	Medical device	Class IIb
Bar code	3661041009291	HS Code	90189084
Warranty (months)	120	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	40	Bare Unit / Width (cm)	11
Packed unit / Height (cm)	15	Bare Unit / Height (cm)	30
Packed unit / Length (cm)	40	Bare Unit / Length (cm)	25.5
Packed unit / Weight	3.155	Bare Unit / Weight (kg)	3.155
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	50
Minimum humidity	1	Maximum humidity	95
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 23:55:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)