

12 GLASSES FOR SPENGLER MANOMETERS VAQUEZ-LAUBRY CLASSIC/LIAN CLASSIC



Ref. 515273



 30/08/2024 19:19:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

12 GLASSES FOR SPENGLER MANOMETERS VAQUEZ-LAUBRY CLASSIC/LIAN CLASSIC



Ref. 515273

Features

CE Marking	NONE	Medical device	No
Bar code	3700446031787	HS Code	901890100000L
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	3
Packed unit / Length (cm)	5	Packed unit / Weight	0.052
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:19:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr