

LIAN NM CLINIC / LARGE ADULT (L)

Ref. 520660



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE

30/08/2024 19:34:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LIAN NM CLINIC / LARGE ADULT (L)



Ref. 520660

Ambidextrous blood pressure monitor, of a robust construction designed for intensive use

Description

High resistance model with Clinic® cuff, quick coupling and clamp for liquid decontamination. Robust, chromium-plated metal casing with foam manometer cover. This model is especially appreciated by emergency services. Supplied in a nylon pouch.

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446024826	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.662
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:34:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr