

MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)



Ref. 522708



CE
0123

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
2
2 ANS
GUARANTEE

30/08/2024 19:40:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)



Ref. 522708

SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.

Description

HIGH MEASUREMENT
 ACCURACY• Double tubing for linear pressure rise

COMFORT• Ambidextrous•
 Remotely readable measurement•
 2.5-meter spiral tubing• Large
 volume ergonomic pear for quick inflation•
 Storage basket behind the dial
 • Comes with cuff wheeled base

• Comes with kit of 3 cuffs (S, M, L), nylon

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3700446031169	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.36
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:40:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr