

AC ADAPTER FOR AUTOTENSIO TENSONIC BPM




Ref. 527184



 30/08/2024 19:44:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

AC ADAPTER FOR AUTOTENSIO TENSONIC BPM



Ref. 527184

Description

AC adapter for electronic blood pressure monitors TENSONIC and AUTOTENSIO arm from SPENGLER

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446022839	HS Code	90330000
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	20
Packed unit / Length (cm)	15	Packed unit / Weight	0.081
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:44:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr