

GREY SPENGLER (M) CUFF (22-36CM) FOR TENSONIC & AUTOTENSIO



Ref. 527463

CE
0197



 30/08/2024 19:45:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

GREY SPENGLER (M) CUFF (22-36CM) FOR TENSONIC & AUTOTENSIO



Ref. 527463

Description

Accessory for ES-60® AUTOTENSIO® TENSONIC® & Effitens®

Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | CE0197 | Medical device | Class I |
| Bar code | 3700446023195 | HS Code | 90262040 |
| Warranty (months) | 6 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.118 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

30/08/2024 19:45:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr