

STATION DE DIAGNOSTIC SPOT-CHECK VS PNI



Ref. 546000



30/08/2024 19:47:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STATION DE DIAGNOSTIC SPOT-CHECK VS PNI



Ref. 546000

Small and light. Easy to transport thanks to its base with wheels.

Description

- Large 10.9 cm screen for good readability•
- Equipped with a long battery life•
- Large memory to record 100 patients
- Three measurement modes: infants, children and adults

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446035822	HS Code	9018199560
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.489
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:47:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr