

STATION DE DIAGNOSTIC SPOT-CHECK 2VS PNI AVEC SOCLE A ROULETTES



Ref. 546002

CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
1
YEAR
GUARANTEE



30/08/2024 19:48:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STATION DE DIAGNOSTIC SPOT-CHECK 2VS PNI AVEC SOCLE A ROULETTES



Ref. 546002

Small and light. Easy to transport thanks to its base with wheels.

Description

- Large 10.9 cm screen for good readability
- Equipped with a long battery life
- Large memory to record 100 patients
- Three measurement modes: infants, children and adults

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446035839	HS Code	9018199560
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	60
Minimum humidity	10	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:48:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr