

SPO2 SENSOR ADULT FOR SPOT-CHECK 2VS



Ref. 546010




CE
0123



 30/08/2024 19:49:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPO2 SENSOR ADULT FOR SPOT-CHECK 2VS



Ref. 546010

Description

SPO2 adult sensor for SPENGLER spot-check station 2VS.

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446029609	HS Code	90181990
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.053
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:49:44 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr