

## SPO2 SENSOR NEONATE FOR SPOT-CHECK 2VS



Ref. 546012



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**II<sub>b</sub>**  
MEDICAL DEVICE

30/08/2024 19:50:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SPO2 SENSOR NEONATE FOR SPOT-CHECK 2VS



Ref. 546012

### Description

SPO2 neonate sensor for SPENGLER spot-check station 2VS.

### Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446029623	HS Code	90181990
Warranty (months)	6	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.07
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 19:50:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)