


SPO2 NEONATE SENSOR FOR 4VS MONITOR



Ref. 551632



CE
0537

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

30/08/2024 20:03:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SPO2 NEONATE SENSOR FOR 4VS MONITOR



Ref. 551632

Description

SPO2 Neonate Sensor for SPENGLER 4 VS Vital Signs Monitor

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031190	HS Code	90181990
Warranty (months)	6	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.086
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:03:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr