

## CONE SINGLE 107 X 14CM (42 INCH) NO PINCH CUFF (THIGH - XL)



Ref. 554058



30/08/2024 20:44:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CONE SINGLE 107 X 14CM (42 INCH) NO PINCH CUFF (THIGH - XL)



Ref. 554058

### Description

Single conical cuffs can be used for a wide spectrum of surgical procedures.

Armbands with zero-pinch technology: the innovative design of the withers of the withers DTS 3000 allows uniform and homogeneous inflation over the entire circumference of the cuff. Better wrapped, the limb requires less pressure to stop blood circulation, reducing here too the risk of complications.

### Features

CE Marking	CE	Medical device	Class I
Bar code	28809315676267	HS Code	9018 90 2090
Warranty (months)	12		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	14
Packed unit / Height (cm)	nc	Packed unit / Length (cm)	nc
Bare Unit / Length (cm)	107	Packed unit / Weight	0.2
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	60
Minimum humidity	5	Maximum humidity	90
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 20:44:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)