

CONE DOUBLE 70 X 10CM (27.5 INCH) NO PINCH CUFF (THIGH - M)



Ref. 554063



30/08/2024 20:45:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CONE DOUBLE 70 X 10CM (27.5 INCH) NO PINCH CUFF (THIGH - M)



Ref. 554063

Description

Double conical cuffs are the best option for long-term surgery to minimize any injury.

Armbands with zero-pinch technology: the innovative design of the withers of the withers DTS 3000 allows uniform and homogeneous inflation over the entire circumference of the cuff. Better wrapped, the limb requires less pressure to stop blood circulation, reducing here too the risk of complications.

Features

CE Marking	CE	Medical device	Class I
Bar code	28809315676472	HS Code	9018 90 2090
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	nc	Packed unit / Length (cm)	nc
Bare Unit / Length (cm)	70	Packed unit / Weight	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	60
Minimum humidity	5	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:45:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr