

LIGHT BULB FOR CLAR MIROR 3.5V (PREVIOUS MODEL 570101) SPENGLER



Ref. 570105



 30/08/2024 20:08:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LIGHT BULB FOR CLAR MIRROR 3.5V (PREVIOUS MODEL 570101) SPENGLER



Ref. 570105

Description

Spare part for Clar mirror

Features

CE Marking	NONE	Medical device	No
Bar code	3700446026936	HS Code	901890800009J
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:08:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr