

## LED BULB FOR CLAR MIRROR AND CHARDON MIRROR




Ref. 570155



 30/08/2024 20:09:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LED BULB FOR CLAR MIRROR AND CHARDON MIRROR



Ref. 570155

### Description

LED bulb for clar mirror and SPENGLER thistle mirror

### Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446000073	HS Code	90189084
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	1	Packed unit / Height (cm)	1
Packed unit / Length (cm)	2	Packed unit / Weight	0.002
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 20:09:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)