

REPLACEMENT HEAD BLUEBERRY FOR OTOSCOPE SMARTLED 5500




Ref. 570518B

 30/08/2024 20:14:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REPLACEMENT HEAD BLUEBERRY FOR OTOSCOPE SMARTLED 5500



Ref. 570518B

Features

CE Marking	NONE	Medical device	No
Bar code	3700446041779	HS Code	90189042
Country	Pakistan		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.046		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:14:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr