

## KERN MGD SCALE, GLASS TRAY, LCD SCREEN, MAX REACH 250 KG

# KERN®

Ref. EKMGD02

GARANTIE  
2 ANS  
2 YEARS  
GUARANTEE



 31/08/2024 00:35:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## KERN MGD SCALE, GLASS TRAY, LCD SCREEN, MAX REACH 250 KG



Ref. **EKMGD02**

### Features

CE Marking	NONE	Medical device	No
Bar code	3361040000199	HS Code	84231090
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

31/08/2024 00:35:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)