

ALUMINIUM BOX 170 X 70 X 20 MM - GREY



Ref. L30200



 30/08/2024 20:38:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ALUMINIUM BOX 170 X 70 X 20 MM - GREY



Ref. L30200

Description

170 x 70 x 20 mm grey

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446000455	HS Code	9018908500
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17
Packed unit / Weight	0.105		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:38:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr