

## THERMOMETRE TEMPO 10 BLUE

Ref. SP1620/



 30/08/2024 20:40:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## THERMOMETRE TEMPO 10 BLUE



Ref. **SP1620/**

### Description

The TEMPO 10 thermometer is a thermometer showing body temperature in 10 seconds•

- Fever alarm•
- Automatic shutdown
- Water resistant mouthpiece•
- Last measurement memory•
- Measurement in 10 seconds
- Comes with a storage case and button cell battery (1.5V)

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446022846	HS Code	90251920
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.051
Bare Unit / Weight (kg)	0.051		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 20:40:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)