

THERMOMETER TEMPO 10 FLEX GREEN



Ref. SP1643/



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

30/08/2024 20:40:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

THERMOMETER TEMPO 10 FLEX GREEN



Ref. **SP1643/**

Description

The TEMPO 10 FLEX thermometer is a thermometer indicating body temperature in 10 seconds and has a flexible mouthpiece•

- Fever alarm•
- Automatic shutdown•
- Water resistant and flexible mouthpiece•
- Last measurement memory•
- Measurement in 10 seconds•
- 2 colors available green and gray
- Comes with a storage case and button battery (1.5V)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446022877	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	6	Packed unit / Height (cm)	2
Packed unit / Length (cm)	15	Packed unit / Weight	0.052
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:40:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr