

## BOX STORE OFFER: 15 BEST-SELLERS AND POS OFFERED



Ref. ZPROMO0407



**LE PACK**  
**15 Produits**  
**1 Présentoir**  
carton à monter  
(L50cm x H125 x P47,2 cm)  
**3 factices\***



31/08/2024 00:32:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## BOX STORE OFFER: 15 BEST-SELLERS AND POS OFFERED



Ref. ZPROMO0407

The best of fixed and mobile DORO phones in a store display

### Description

We have selected for you 8 references, all Best-sellers (current / and in the making)  
 A Pack that includes:>> 15 Products +  
 1 Cardboard display to mount  
 + 3 dummies\*

In detail, the products contained in the Offer: 3 2G compact mobiles: doro

- 1380 - Black or Red
- 3 Compact 4G mobiles: doro 5860, Graphite or Black & White
- 2 2G clamshell mobiles: doro 2424 - Graphite & white or Red & white 3
- 4G clamshell mobiles: doro 6880 -Black & White or Red & White 1 Cordless landline phone doro Phone Easy 100w

- White Large Display 30db 1
- Cordless landline phone DUO doro phone easy 100w, white 1 Wired Phone doro Phone Easy 331ph 3 memo-photo keys
- 1 Wired Phone doro Magna 4000, 60dB: the highest amplification in the doro

range

### Features

CE Marking	NONE	Medical device	No
Bar code	3661040221243	HS Code	85171800
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

31/08/2024 00:32:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
 13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)