

## OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY



Ref. 221052



30/08/2024 20:55:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY



Ref. 221052

For SpO2 spot check.

### Description

#### HIGH MEASUREMENT

ACCURACY• Plethysmography curve display• High sensitivity sensor

#### EASY TO USE AND

CONVENIENT• Compact and lightweight• Switch on at the touch of the clock• Automatic shutdown after 5 seconds of non-use• Automatic display rotation• Simultaneous display in less than 10 sec. :  
 - pulsed oxygen saturation (% SpO2),  
 - heart rate,  
 - the intensity of the pulsatile signal (graph bar),  
 - the plethysmography curve.  
 • TFT display• Low battery indicator• Customizable display: 6 display modes.  
 • 5 levels of brightness

#### ADDITIONAL

INFORMATION• Comes with a wrist strap, a cover and 2 AAA LR03 batteries

### Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446041809	HS Code	90181910
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.135
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-40	Max temperature	60
Minimum humidity	5	Maximum humidity	95
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 20:55:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
 13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)