

PULSE OXYMETER OXYGO GREY



Ref. 221350/



30/08/2024 20:56:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PULSE OXYMETER OXYGO GREY



Ref. 221350/

For adults and children

Description

Suitable for healthcare professionals wishing to have an oximeter that is easy to use and quick to implement, for use in the office or on the go.

VERY HIGH MEASUREMENT
 ACCURACY• Child-friendly technology•
 High sensitivity sensor

EASY TO USE AND
 CONVENIENT• Measurement in 10 seconds•
 Shock resistant
 • Very long autonomy•
 Desaturation indicator•
 Very easy to clean because developed without any roughness

ADDITIONAL INFORMATION
 • Comes with a wrist strap, a cover and 1 AAA LR03 battery

Features

CE Marking	CE0197	Medical device	Class IIa
Bar code	3700446021719	HS Code	90181910
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.149
Bare Unit / Weight (kg)	0.04		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:56:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr