

SINGLE FULL SINGLE AUTOCLAVABLE CUFF 52 X 7CM (20.8 INCH)



Ref. 554075



30/08/2024 20:49:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SINGLE FULL SINGLE AUTOCLAVABLE CUFF 52 X 7CM (20.8 INCH)



Ref. 554075

Description

This cuff for the withers DTS 3000 is completely autoclavable and consists of a special fabric coating and suitable Velcro, an internal silicone bladder and silicone hoses that can be fully sterilized in an automatic autoclave machine.

Armbands with zero-pinch technology: the innovative design of the withers of the withers DTS 3000 allows uniform and homogeneous inflation over the entire circumference of the cuff. Better wrapped, the limb requires less pressure to stop blood circulation, reducing here too the risk of complications.

Features

CE Marking	CE	Medical device	Class I
Bar code	28809315676601	HS Code	9018 90 2090
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Packed unit / Length (cm)	nc
Bare Unit / Length (cm)	52	Packed unit / Weight	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:49:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr