

## DTS-3000 STAND, 102CM X 54CM

Ref. 554115



30/08/2024 20:51:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DTS-3000 STAND, 102CM X 54CM



Ref. 554115

### Features

CE Marking	CE	Medical device	No
Bar code	18809315676182	HS Code	9018 90 2090
Warranty (months)	12		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	56	Packed unit / Height (cm)	23
Packed unit / Length (cm)	56	Packed unit / Weight	11.6
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 20:51:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)