

LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 WITH CHARGING DOCK



Ref. 579530



30/08/2024 20:21:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LED FIBER OPTIC OTOSCOPE - SMARTLED 5500 WITH CHARGING DOCK



Ref. 579530

Rechargeable LED fiber optics with Smartled 5500-R® charger stand

Description

- Cutting-edge LED technology:- A powerful luminance of 50K lumens,- A colorimetry of 4000K, the perfect range for observation, - A CRI (Color Rendering Index) of 95 out of 100 for perfect color reproduction.
- Optical fiber for direct, homogeneous and glare-free illumination of the ear canal and eardrum• Magnifying magnifying glass x3

EASY TO USE AND PRACTICAL•

- Ergonomic handle•
- ABS and stainless steel head for secure fixation•
- Clip with switch•
- Crimped bulb for greater reliability•
- LED bulb with a lifespan of 55000 h

ADDITIONAL INFORMATION

- Comes with rigid housing consisting of 2 LR6 AA batteries, 8 speculums (2.5mm and 4mm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446032937	HS Code	94029010
Warranty (months)	36	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	12
Packed unit / Length (cm)	23	Packed unit / Weight	0.766
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:21:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr