

DISPOSABLE OXYGEN MASK WITH TUBING HIGH CONCENTRATION PEDIATRIC



Ref. 580873


CE
0120



 30/08/2024 20:24:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DISPOSABLE OXYGEN MASK WITH TUBING HIGH CONCENTRATION PEDIATRIC



Ref. **580873**

Description

Single-use medium or high concentration administration mask.
High concentration: the mask can deliver oxygen concentrations close to 100% for flow rates of 10-15 l/min. With adult or child tubing.

Features

CE Marking	CE0120	Medical device	Class IIa
Bar code	5030267051779	HS Code	90192000
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Packed unit / Height (cm)	6
Packed unit / Length (cm)	15	Packed unit / Weight	0.068
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:24:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr