

SELF PIERCING PEN

Ref. 745011



30/08/2024 20:30:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SELF PIERCING PEN



Ref. 745011

Description

- Lancing pen for use with SPENGLER lancets.
- Packaging: Unit

The lancing pen must be used with compatible sterile lancets to extract a blood sample taken from the finger, palm or forearm to test blood glucose or other tests using small amounts of blood.

Find our complete kit Glucometer SPENGLER TB 100 (ref. 745000)

LPP code for reimbursement of the lancing pen: 611139672

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041892	HS Code	90183190
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	2	Bare Unit / Width (cm)	0.8
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	3.5
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	12
Packed unit / Weight	0.22	Bare Unit / Weight (kg)	0.02
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Code LPP	611139672	Packed unit / Quantity	1

30/08/2024 20:30:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr