

VINYL RESUSCITATORS CHILDREN SINGLE USE (MASK 3)



Ref. 810104



CE
0120



 30/08/2024 20:32:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

VINYL RESUSCITATORS CHILDREN SINGLE USE (MASK 3)



Ref. 810104

Description

Emergencies - resuscitation - first aid. Flexible, lightweight and ergonomic, these resuscitators maintain safe and effective ventilation, as much in cases requiring resuscitation as in cases where breathing is low. These systems can be connected to an oxygen supply.

Features

CE Marking	CE0120	Medical device	Class IIa
Bar code	5030267073245	HS Code	90192000
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.55
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:32:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr