

VINYL MASK SINGLE USE, CHILDREN, WHITE, T2



Ref. 860131




CE
0120



 30/08/2024 20:34:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

VINYL MASK SINGLE USE, CHILDREN, WHITE, T2



Ref. 860131

Description

Emergencies - resuscitation - first aid. Can be fitted to the manual insufflators. These facial resuscitation masks ensure administration of oxygen and anesthetic vapours from anesthesia administration devices. They can be used in spontaneously breathing patients or to ensure controlled ventilation. Their specially designed shape means they fit all morphologies.

Features

CE Marking	CE0120	Medical device	Class IIa
Bar code	5030267078912	HS Code	90192000
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.025
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:34:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr