

## COMBAT APPLICATION TOURNIQUET, BLACK COLOR

Ref. STGARR01



31/08/2024 00:10:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## COMBAT APPLICATION TOURNIQUET, BLACK COLOR



Ref. **STGARR01**

Tourniquet very used by the military and rescuers.

### Description

- Strong and resistant
- Ready to use to intervene quickly in case of hemorrhage.
- Easy to apply, it can be used with one hand, and it allows rescuers to quickly free themselves to go see other victims.
- Suitable for arms and leg (including thigh)
- Dimensions: 90 x 3.8 cm
- A label is present to write the time of the withers inserted.
- Black color, it also exists in Orange

### Features

CE Marking	CE	Medical device	Class I
Bar code	3361040000205	HS Code	3006100000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.078
Bare Unit / Weight (kg)	0.078		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

31/08/2024 00:10:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)