

## FETAL AND VASCULAR DOPPLER WITHOUT PROBE



Ref. 743301



16/09/2024 10:27:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## FETAL AND VASCULAR DOPPLER WITHOUT PROBE



Ref. 743301

**Small and lightweight, this pocket Doppler is the essential tool for medical practices and communities**

### Description

For prenatal ultrasound and examination of blood vessels. Choose the probes suited to your requirements. Large colour LCD screen displaying hear t rate, battery indicator, pulse wave intensity and probe frequency. Supplied with two LR6 AA batteries. Ergonomic angled probes for a comfortable use.

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446004316	HS Code	90181990
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.43
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	55
Minimum humidity	5	Maximum humidity	93
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

16/09/2024 10:27:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)