

GUEDEL CANNULA SINGLE USE STERILE ISO 9.0 ORANGE, N°3)




Ref. 580842



 16/09/2024 10:10:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

GUEDEL CANNULA SINGLE USE STERILE ISO 9.0 ORANGE, N°3)



Ref. 580842

Description

Single-use cannula in individual sterile bag.
The Guedel cannula is used to keep the airway of an unconscious patient open.

Features

CE Marking	NONE	Medical device	Class I
Bar code	5032670506808	HS Code	90183900
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.013
Bare Unit / Weight (kg)	0.013		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

16/09/2024 10:10:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr