

TONGUE PULLER, ORANGE

Ref. 580864



16/09/2024 10:13:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

TONGUE PULLER, ORANGE



Ref. 580864

Description

Made of autoclavable composite material at 134°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040212739	HS Code	90200000
Warranty (months)	60	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.01
Bare Unit / Weight (kg)	0.01		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

16/09/2024 10:13:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr