

LUXASCOPE DERMATOSCOPE LED 3.7 V (RECHARGEABLE), INCL. USB-RECHARGER

LUXAMED®

Ref. C1.526.114



Rechargeable USB



17/09/2024 17:48:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LUXASCOPE DERMATOSCOPE LED 3.7 V (RECHARGEABLE), INCL. USB-RECHARGER

LUXAMED®

Ref. **C1.526.114**

Innovative dermatoscope in pocket size, rechargeable USB.

WITH scale.

Description

- Illumination: LED-Ring with 6 LEDs in concentric structure for a highly detailed image, homogenous light
- Maintenance-free, because no lamp replacement necessary
- Removable contact plate for easy cleaning
- Stepless zoom area with 10 x magnification for a distortion-free image
- Optimized design of head
- Can be dimmed
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours by constant current drive
- Intensity of illumination: 10.000 Lux
- Colour temperature: approx. 4.200 K
- Two-component handle (aluminium / fibre-glass reinforced), also combinable with LuxaScope LED Otoscope- and Ophthalmoscope head
- Electronic ON/OFF push-button-switch
- Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive
- Incl. recharger with USB plug (adapter EU/UK/US)

Features

CE Marking	NONE	Medical device	No
Part code	04260381963448	Warranty (months)	24

DIMENSIONS AND WEIGHT

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr

OTHER FEATURES

Colors

Black