

DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

LUXAMED®

Ref. D1.312.912


CE



 17/09/2024 20:24:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

LUXAMED®

Ref. **D1.312.912**

Penlight for eye and throat diagnostic., long life technology.

Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963349	Warranty (months)	24
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Kids		

📅 17/09/2024 20:24:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr