



Ref. 507219



17/09/2024 16:25:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

Ref. 507219

Spengler takes care of newborns!

Description

The dual Magister Neonatal est dedicated to newborn auscultation.

With its dual head chestpiece equipped of a floating diaphragm on the large side, Magister Neonatal is perfectly adapted to babies for listening high and low frequencies.

Its PVC non-chill rings on both sides of the chestpiece allows a perfect comfort for the newborns.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446042325	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Bare Unit / Length (cm)	76
Packed unit / Weight	0.25		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

 17/09/2024 16:25:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr