

TEMPO MAX HYPOTHERMIC DIGITAL THERMOMETER SPENGLER

Ref. THE002



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 ANS
2 YEARS
GUARANTEE

17/09/2024 15:49:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

TEMPO MAX HYPOTHERMIC DIGITAL THERMOMETER SPENGLER

Ref. THE002

=====

Description

=====

Features


CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392219271	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.3
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.02
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES

 17/09/2024 15:49:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr