

BAG OF 10 UNITS ECG ADAPTER PRESSURE BUTTON



Ref. 260023



 01/10/2024 20:30:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BAG OF 10 UNITS ECG ADAPTER PRESSURE BUTTON



Ref. 260023

=====

Description

Kit of 10 adapters for snap ECG electrodes CARDIOMATE

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030704	HS Code	90330000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.06
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 01/10/2024 20:30:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr