

STAINLESS STEEL STIRRUPS FOR VOG ELECTRIC EXAMINATION COUCH


Ref. OV0013



03/10/2024 21:58:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

STAINLESS STEEL STIRRUPS FOR VOG ELECTRIC EXAMINATION COUCH



Ref. OV0013

Discover our stainless steel stirrups for electric sofa

Features

CE Marking	NONE	Medical device	No
Bar code	3661040094342	HS Code	94029000
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/10/2024 21:58:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr