

STATION DE DIAGNOSTIC SPOT-CHECK 2VS PNI+SPO2



Ref. 546001



09/10/2024 00:03:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STATION DE DIAGNOSTIC SPOT-CHECK 2VS PNI+SPO2



Ref. 546001

Small and light. Easy to transport thanks to its base with wheels.

Description

- Large 10.9 cm screen for good readability•
- Equipped with a long battery life•
- Large memory to record 100 patients
- Three measurement modes: infants, children and adults. Oxygen saturation module.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6941900603387	HS Code	9018199560
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.439
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	60
Minimum humidity	10	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		

09/10/2024 00:03:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr